



Improvement in Bag-Mask Ventilation Performance After Training With a Novel Terminal Feedback Manikin System

Angelo A. Salvucci, MD, Benjamin T. Squire, MD, Amy H. Kaji, MD, PhD, James T. Niemann, MD, David G. Chase, MD
Ventura County Medical Center, Harbor-UCLA Medical Center



BACKGROUND

- Bag-mask ventilation (BMV) is commonly used in emergency and prehospital care.
- With the de-emphasis of endotracheal intubation, BMV has become increasingly important in cardiac arrest and prehospital management of traumatic brain injury.
- BMV is considered a basic skill, often taught only by demonstration and brief practice.
- Poorly performed BMV is common, resulting in brief, high-volume inspirations delivered at a rapid rate.
- Consequences include gastric inflation, an increase in intrathoracic pressure with resultant hypotension, and hypocarbia reducing cerebral blood flow.
- Despite poor performance and known complications, training methods have remained largely unchanged.
- A well known and effective method to learn motor skills is terminal feedback - performance assessment provided after task completion.

OBJECTIVE

- To measure adequacy of performance in bag-mask ventilation in each of three parameters: duration of inspiration (1 second), volume of inspiration (500-600 mL), and duration of breath-to-breath interval (5-6 seconds).
- To evaluate the performance improvement impact of visual feedback training.



Ambu SmartMan® Feedback Display

METHODS

Design: Prospective pre-post interventional study.
Subjects: Twenty-eight volunteer on-duty EMTs and paramedics performed bag-mask ventilation on a SmartMan® adult CPR manikin system from Ambu, Inc. and remained blinded to the results.

- Subjects were then trained for 1 minute using the SmartMan continuous terminal feedback (bandwidth knowledge of results) training system, utilizing a color bar display to indicate specified ranges for “acceptable”, “too high” and “too low” on each of the three parameters. A second blinded assessment was then done for all subjects.
- To account for potential correlations among ventilations for each subject, a cluster analysis was performed using generalized estimating equations and an exchangeable covariance matrix; proportional odds ratios (OR) with 95% confidence intervals were calculated.

Outcome measures:

- Adequacy of duration of inspiration, volume of inspiration, and duration of breath-to-breath interval.
- Performance differences before and after training intervention.




SmartMan Training/Testing

RESULTS

- Twenty-eight subjects performed a total of 672 total breaths, 336 per trial.
- Individual measurements were considered correct if performed to within 20% of the AHA Guidelines for ventilations.
- Acceptable ranges: Duration of inspiration 0.8-1.2 seconds, volume of inspiration 400-700 mL, breath-to-breath interval 4-7 seconds.
- Correctly performed duration of inspiration increased from 34% of ventilations to 80%.
- Proportion delivered too rapidly decreased from 50% to 15%.
- Correct volume of inspiration increased from 63% to 96%.
- Correct breath-to-breath interval increased from 63% to 93%.
- All pre-post intervention differences were significant (P<0.0001).

BMV Performance (% correct, 95% CI)



CONCLUSIONS

- There is accumulating evidence that deficiencies in BMV skills worsen patient outcomes.
- Initial BMV was performed inconsistently over a wide range of durations and volumes.
- Too-rapid inspiration was the most common error, occurring in nearly 50% of breaths.
- A brief visual terminal feedback training session substantially improved performance, with a large majority of breaths delivered within acceptable ranges.
- Correct inspiratory rate and volume with bag-mask ventilation appears to be a learnable skill.
- Terminal feedback methods of training can optimize motor skills learning and performance.

REFERENCES

- Aufderheide TP, Lurie KG. Death by hyperventilation: a common and life-threatening problem during cardiopulmonary resuscitation. *Crit Care Med* 2004;32 (Suppl):S345-S351.
- Pepe PE, Roppolo, LP: The detrimental effects of ventilation during low-blood-flow states. *Curr Opin Crit Care*. 2005 Jun;11(3):212-8.
- Davis DP, Idris AH, Sise MJ, et al: Early ventilation and outcome in patients with moderate to severe traumatic brain injury *Crit Care Med*. 2006 Apr;34(4):1202-8.
- Warner KJ, Cuschieri J, Copass MK, Jurkovich GJ, Bulger EM The impact of prehospital ventilation on outcome after severe traumatic brain injury. *J Trauma*. 2007 Jun;62(6):1330-8.
- Aufderheide TP, Sigurdsson G, Pirralo RG, et al. Hyperventilation-induced hypotension during CPR. *Circulation* 2004; 109:1960-1965.
- Spooner BB, Fallaha JF, Kocierz L: An evaluation of objective feedback in basic life support (BLS) training. *Resuscitation*. 2007 Jun;73(3):417-24.